



## HIGH LEVEL DISINFECTANT

# ANIOXYDE 1000 LD and BANDELETTES ANIOXYDE 1000 LD (Test Strips)

### HIGH LEVEL DISINFECTANT FOR MEDICAL DEVICES

- Shelf life before use: 24 months
- Patented formula
- High level disinfectant from PHERA® system, active against spores of bacteria
- Extemporaneous production of peracetic acid
- Alternative to glutaraldehyde
- Aldehyde free: no fixation of proteins
- Possible control of bath efficacy with test-strips




# ANIOXYDE 1000 LD

High level disinfectant for medical devices


## PRODUCT CHARACTERISTICS

- Bath colour: pink after activation
- Absence of acetic acid vapour: better user tolerancy
- Easy use: activator integrated to the 1 L & 5 L can
- Efficacy proved at 1000ppm of peracetic acid
- Use duration of the soaking bath: 14 days maximum (7 days for France) check the conformity of the solution with test-strips
- Wide compatibility with materials, tanks and working table
- Can be used in the presence of a vapour suction system
- The disinfecting solution contains from 2000 ppm of peracetic acid (after activation) to 1000 ppm (threshold limit concentration detected by test-strips)


## INSTRUCTIONS FOR USE




**Pre-disinfection step:** Clean the medical device with a detergent or a pre-disinfectant product. Rinse thoroughly. For endoscopic equipment: brush channels.



**Disinfection step:** Pour the activated solution (see activation procedure on the label/protocol). Cover the basin. Respect the minimum contact time. For endoscopic equipment: make solution circulate inside channels. Do not add any product into Anioxyde 1000 LD.



Rinse thoroughly with sterile or filtered water (0.2 µm) to eliminate any trace of the product. For endoscopic equipment: rinse the internal and external parts of the medical device.



Dry with a single-use towel. Keep the medical device as aseptically as possible until its next use.

## MICROBIOLOGICAL PROPERTIES

CLAIM/CONDITIONS	CONTACT TIME
Bactericidal - According to EN 13727/EN 14561, clean conditions	5 min
Yeasticidal - According to EN 13624/EN 14562, clean conditions	5 min
Mycobactericidal - According to EN 14348/EN 14563, clean conditions	5 min
Fungicidal - According to EN 13624/EN 14562, clean conditions	5 min
Virucidal - According to EN 14476/EN 17111, clean conditions	5 min
Adenovirus - According to EN 14476/EN 17111, clean conditions	5 min
Norovirus - According to EN 14476/EN 17111, clean conditions	5 min
Sporicidal* against Clostridium difficile and Bacillus subtilis according to EN 17126, clean conditions	5 min

## PACKAGING AND ORDERING INFORMATION

PRODUCT	PACKAGING UNIT	ORDER CODE
Anioxyde 1000 LD	4 x 5 L canister + integrated activator	XXXXXX
Anioxyde 1000 LD	12 x 1 L canister + integrated activator	XXXXXX
Anioxyde 1000 LD	Test strips of Peracetic Acid rate 1000 ppm (50 per box)	XXXXXX



## APPLICATION AREA

High level disinfectant for medical devices: surgical and medical equipment, endoscopy equipment and heat-sensitive equipment.

## INSTRUCTIONS FOR USE

Two-component-solution. Ready to use after activation.

**For 5 L canister:** Pour all of the activator (bottle containing 39.1 ml) into the generator (5 L container) and then pour the solution produced into the soaking bath.

**For 1 L canister:** Pour all of the activator (bottle containing 8.3 ml) into the generator (1 L container) and then pour the solution produced into the soaking bath.

Before the first use, wait at least 30 minutes (time required for activation of the product). Carefully clean and rinse the medical device before soaking it in ANIOXYDE 1000 LD (activated solution).

Check the level of peracetic acid with a ANIOXYDE 1000 LD test strip at the beginning of the session, when resuming activity after an interruption and every 4 hours in case of continuous activity. Maximum shelf life of the soaking bath: 14 days (except for France 7 days maximum according to French regulations).

1. Pre-disinfection step: Clean the medical device with a detergent or a pre-disinfectant product. For endoscopic equipment: brush channels.
2. Disinfection step: Pour in the activated solution (see activation procedure on the label/protocol). Keep the soaking bath covered during the contact time. Respect the minimum contact time. For endoscopic equipment: make solution circulate inside channels. Do not add any product into Anioxyde 1000 LD.
3. Rinse thoroughly with sterile or filtered water (0.2 µm) to eliminate all traces of the product. For endoscopic equipment: rinse the internal and external parts of the medical device.
4. Dry with a single-use towel. Keep the medical device as aseptically as possible until its next use.

## INGREDIENTS

Extemporaneous production of peracetic acid from the acetylcaprolactam (PHERA system), Anioxyde 1000 LD generator containing 3% hydrogen peroxide and stabilising agents.

## STORAGE

Dangerous - respect the precautions for use (Drawn up according to the European rules in force regarding the classification and labelling of chemical products. Storage: between +5°C and +25°C. STORE UPRIGHT.

## CERTIFICATES

Medical device. Full antimicrobial activity for Anioxyde 1000 LD is included in the scientific dossier, available on request.

## MEDICAL DEVICE USE SAFELY

Please read label and product information.



FOR PROFESSIONAL USE ONLY.

# BANDELETTES ANIOXYDE 1000 LD (Test Strips)

## INDICATIONS

50 test strips for the level of peracetic acid in the ANIOXYDE 1000 LD reconstituted solution.

## INSTRUCTIONS FOR USE

- Dip the strip for one second, remove it while shaking it to remove excess liquid
- Reading after 10 seconds:  
Uniform black across the entire reaction range - COMPLIANT  
Non-uniform colouring, evidence of white spots, colour soon fades - NON COMPLIANT
- Never interpret a result more than 15 seconds after dipping
- For the frequency of use, refer to the instructions on the ANIOXYDE 1000 LD label

## STORAGE CONDITIONS

Close the container after removing each testing strip. Keep away from moisture. Respect the expiry date shown on the label. Store up to 90 days after opening. Store between +5°C and +25°C.

## PRECAUTIONS FOR USE

Remove the product and its container by following all the warnings of use. Reserved for professional use.




## INSTRUMENT REPROCESSING

### READING AFTER 10 SECONDS:

**Compliant**  
Uniform black across the entire reaction range



**Non compliant**  
Non-uniform colouring, evidence of white spots, colour soon fades



FOR PROFESSIONAL USE ONLY.



**LABORATOIRES ANIOS**

1, Rue de l'Espoir  
59260 Lezennes  
France

Distributed by:

**ECOLAB DEUTSCHLAND GmbH**

Ecolab-Allee 1  
40789 Monheim am Rhein  
Germany



[www.ecolab.com](http://www.ecolab.com)