



AseptoCool

Powerful alkaline cleaner for cold or hot cleaning of milking equipment and milk tanks.

- ▲ Hot or cold cleaning
- ▲ Hard or soft water

APPLICATION

Pre-clean with water. Clean with **AseptoCool** 1.0% (at room temperature, 10 min.) or 0.7% (at 70°C; 10 min.). Rinse with water of drinking quality. To be used 12 to 13 times a week - in rotation with, once or twice a week, **HorolithCool**.



FEATURES AND BENEFITS

AseptoCool is an alkaline chlorinated product with boosted cleaning properties. It is suitable for all kinds of standard milking equipment, including cooling tanks.

The **Cool** products can be used at a wide temperature range including cold washes and in all water hardnesses.

SPECIFICATIONS

AseptoCool	
Physical Form	Clear, light yellow liquid
Active Substance	Sodium hydroxide, Sodium hypochlorite
pH at 20°C:	13 - 14
Density at 20°C (g/cm ³):	1.26 - 1.29
Storage temperature:	-10°C to +40°C

LOGISTICS

Packaging	13 kg	25 kg	265 kg
Product code	2321680	2321690	2321700
EAN code	4028162321685	4028162321692	4028162321708

A Material Safety Data Sheet is available on request for professional users.
A Product Data Sheet is available on request.

Ecolab Europe GmbH
Richtistrasse 7
CH-8304 Wallisellen
www.ecolab.com