



AseptoForce



Powerful chlorinated alkaline cleaner and disinfectant for milking equipment and milk tanks.

- ▲ Strong disinfection
- ▲ Ideal alkalinity due to potassium and sodium hydroxide
- ▲ For challenging conditions

APPLICATION

Pre-clean with water. Clean and disinfect with **AseptoForce** at 1 % (minimum 10 min. at 40 - 50°C). Rinse with water of drinking quality. To be used in daily rotation with **HorolithForce**.



FEATURES AND BENEFITS

AseptoForce is a well-balanced alkaline chlorinated product, containing sodium and potassium hydroxide, with excellent cleaning and disinfection properties. It is suitable for all kinds of standard milking equipment, including cooling tanks. **AseptoForce** is in compliance with the Biocide Product Regulation (No. 528/2012).

The **Force** products are specially developed for challenging cleaning conditions. The combination of dedicated ingredients assures best hygiene results.

SPECIFICATIONS

AseptoForce	
Physical form	Light yellow liquid
Active substance	Sodium hypochlorite
pH at 20°C:	13.0 - 14.0
Density at 20°C (g/cm ³):	1.14 - 1.18
Storage temperature:	-15°C to +30°C

LOGISTICS

Packaging	12 kg	23 kg	70 kg	240 kg
Product code	2331540	2331550	2334640	2331570
EAN code	4028162331547	4028162331554	4028162334647	4028162331578

A Material Safety Data Sheet is available on request for professional users.
A Product Data Sheet is available on request.

Ecolab Europe GmbH
Richtstrasse 7
CH-8304 Wallisellen
www.ecolab.com