



# AseptoRobot Des



Powerful chlorinated alkaline cleaner and disinfectant for automatic milking systems (AMS), standard milking equipment and milk bulk tanks.

- ▲ Designed for all kinds of AMS
- ▲ Disinfectant chlorinated alkaline
- ▲ Designed for warm water cleaning (40 - 70°C)

## APPLICATION

1- AMS: Follow the AMS cleaning routine. To be used aligning with robot cleaning program in rotation with **HorolithRobot**.  
 2- Milking equipment (non AMS): Pre-clean with water. Clean and disinfect with **AseptoRobot Des** at 0.5 - 0.7% (minimum 10 min. at 40°C - 50°C). Rinse with water of drinking quality. To be used in daily rotation with **HorolithRobot**.



## FEATURES AND BENEFITS

**AseptoRobot Des** is an alkaline chlorinated product to clean and disinfect the AMS. It is suitable for all AMS, that do not exceed 70°C during the cleaning program. It can also be used for cleaning and disinfection of other milking equipment, including cooling tanks. **AseptoRobot Des** is in compliance with the Biocide Product Regulation (No. 528/2012).

The **Robot** products are specially developed for use in AMS.

## SPECIFICATIONS

AseptoRobot Des	
Physical form	Light yellow liquid
Active substance	Sodium hypochlorite
pH at 20°C:	13.5 - 14.0
Density at 20°C (g/cm <sup>3</sup> ):	1.18 - 1.22
Storage temperature:	-15°C to +30°C

## LOGISTICS

Packaging	24 kg	245 kg
Product code	2331440	2331450
EAN code	4028162331448	4028162331455

A Material Safety Data Sheet is available on request for professional users.  
 A Product Data Sheet is available on request.

**Ecolab Europe GmbH**  
 Richtstrasse 7  
 CH-8304 Wallisellen  
[www.ecolab.com](http://www.ecolab.com)