

Phytoshield® D

Natural and dynamic Shield to disinfect and protect teats after milking

- ▲ Ultimate disinfecting power of lactic acid
- ▲ Active polymerising Shield
- ▲ Astringent properties accelerating the closing of the sphincter in the teat canal



APPLICATION

Phytoshield® D is used neat. Fill the reservoir of the dipping cup and apply on the teats after milking. Do not wipe. Prevent cows from coming into contact with bedding material until product is dry.



PREMIUM FEATURES

Phytoshield® D contains lactic acid, a notified active substance for disinfecting the skin of teats for veterinary hygiene. This substance contains plant extracts providing immense disinfection and hypoallergenic properties.



The polymerizing Shield formed by Phytoshield® D protects the teats from physical and microbial soils from one milking session to the next. The formulation enables easier cleaning at the next milking.

BENEFITS

Phytoshield® D provides dual protection at the entrance of the teat canal with its superior shield technology and the action of lactic acid. The polymers in Phytoshield® D have been specially selected to let the skin breathe without blocking metabolic exchanges. Thanks to its unique stringent properties, Phytoshield® D accelerates the closure of the teat canal reinforcing its action against any difficult environment.

LOGISTICS

Packaging:	10kg	20kg	60kg	210kg
Product code:	2289110	0000000	0000000	0000000
EAN code:	402816 0000000	402816 0000000	402816 0000000	402816 0000000
Palletization:	60 cans	24 cans	6 units	4 units



SPECIFICATIONS

Product name :	Phytoshield® D
Physical Form	Liquid
Active Substance	Lactic acid
Fragrance	Mint
pH Value	3,75 ± 0,15
Storage Temp.	0°C – 30°C